

# Filter für Chemikalien-Laborabzüge ChemFAST Elite

Alle Filter werden aus hochwertiger Kokosnussschalen- oder Steinkohle hergestellt. Die amorphe Struktur und hohe Absorptionseffizienz der Aktivkohle aus Kokosnussschalen wird durch die kontrollierte Hitzeoxidation erreicht.

Typ	Anwendung Filter	Für	Menge pro VE	Bestell Nr.
A/C	Lösungsmittel	ChemFAST Elite 06/12/15	1	4.691 992
A/C	Lösungsmittel	ChemFAST Elite 09/15/18	1	4.692 001
SULF	Säuren, Sulfide RH >85 %	ChemFAST Elite 06/12/15	1	4.691 993
SULF	Säuren, Sulfide RH >85 %	ChemFAST Elite 09/15/18	1	4.692 002
UR	Säuren, Sulfide RH <85 %	ChemFAST Elite 06/12/15	1	4.691 994
UR	Säuren, Sulfide RH <85 %	ChemFAST Elite 09/15/18	1	4.692 003
AM	Ammoniak	ChemFAST Elite 06/12/15	1	4.691 995
AM	Ammoniak	ChemFAST Elite 09/15/18	1	4.692 004
CYAN	Blausäure	ChemFAST Elite 06/12/15	1	4.691 996
CYAN	Blausäure	ChemFAST Elite 09/15/18	1	4.692 005
MER	Quecksilber	ChemFAST Elite 06/12/15	1	4.691 997
MER	Quecksilber	ChemFAST Elite 09/15/18	1	4.692 006
ACR	radioaktives Jod	ChemFAST Elite 06/12/15	1	4.691 998
ACR	radioaktives Jod	ChemFAST Elite 09/15/18	1	4.692 007
FOR	Formaldehyd	ChemFAST Elite 06/12/15	1	4.691 999
FOR	Formaldehyd	ChemFAST Elite 09/15/18	1	4.692 008